

CENTRINO PUMP FOR FILTRATION

Product Code: EEFP1H

HSN Code: 90318000

The Centrino Filtration Pump is a highly efficient and it saves 70% Desk space, The compact design of the filtration system can save much Desk space by simplifying the need for a separate manifold, pump, and excess tubing. Improve operation time, direct drainage design saves a lot of time to empty waste bottle and make the filtration more efficient.

FEATURES:

- ☑ Dual Power Operated, Battery Back 1 hr 30mins. And DC Power.
- ☑ Data Easy to sanitize with SS316L Body
- ☑ Replaceable funnel adaptor design makes Centrino Pump compatible with various kinds of filter holder User.
- ☑ Improve the Operation time.
- ☑ Compact design saves 70%Desk space.



TECHNICAL PARAMETERS

Model	Centrino Filtration Pump
Max. liquid flow rate	260ML Per min
Power	DC 24 V, Built-in Lithium battery
Motor	DC 24 V Pump
Hose barb	¼ Inch (6 mm)
Material contacts with	SS316L, Silicon Seal, PP Filtrate
Dimension	125mm(L)*160mm(W)*150mm(H)
Weight	1.6Kg
Number of Branch	Single Head
Filter holder	MF47 Funnel, PSU Funnel