

HEATING BLOCK INCUBATOR

Product Code : EEDI60W

HSN Code : 90181290

TECHNICAL PARAMETERS

Product Name	TAL-M2 Dry Heat Incubator (Double Modules)
Temperature Range	Ambient +5°C ~ 150°C
Time Required	≤ 30 min (from 20°C rise to 150°C)
Temperature Stability @100~150	≤ ±1°C
Temperature Stability @40~100	≤ ±0.5°C
Temperature Uniformity 40	0.3°C
Displayed Precision	0.1°C
Temperature Control time	99h59min
Highest Temperature	150°C
Power	400W
Dimension	260x290x95mm 360x220x95mm

