

PURE STEAM QUALITY TEST KIT

Product Code : EPL STQT

HSN Code : 84219900

Elena Electronics Pacific Pvt. Ltd Design to deliver the best products to our customers. To achieve this and try to get the experts opinion for the world class products For Steam Quality Test Kit. Pure Steam Quality plays a vigorous role in steam sterilization of load in Pharmaceutical and Biological Industry. International regulatory standards such as EN 285, EN 17665 & HTM2010 and ISO 111343 recommends steam quality testing and has defined parameters for compliance of clean steam. It is defined as the steam its condensate meets all parameters of WFI. It must meet the basic following physical properties.

- Non-Condensable gases < 3.5%
- Dryness Value >0.9%
- Super Heat Value <25oC

Elena offers Pure Steam Quality Test Kit that enables engineer to provide dependable and wellorganized testing without needs of bulky equipment.

BEST FEATURES

- Designed for quick setup
- Easy to use
- Efficient
- Self-contained
- Comprehensive
- Suitable



Elena Electronics Pacific Pvt. Ltd, Pure Steam Quality Test Kit help you identify steam conditions & meeting regulatory requirements quickly & easily. Complete solution with Validation, Training & Documentation.