

## CENTRINO PUMP FOR FILTRATION – 3 HEAD

Product Code: EEFP3H

HSN Code: 90318000

The Centrino Filtration Pump is a highly efficient and it saves 70% Desk space, The compact design of the filtration system can save much Desk space by simplifying the need for a separate manifold, pump, and excess tubing. Improve operation time, direct drainage design saves a lot of time to empty waste bottle and make the filtration more efficient.

### FEATURES:

- Dual Power Operated, Battery Back 1 hr 30 mins. And DC Power.
- Data Easy to sanitize with SS316L Body
- Replaceable funnel adaptor design makes Centrino Pump compatible with various kinds of filter holder User.
- One - Year Warranty.
- Improve the Operation time.
- Compact design saves 70%Desk space.



### SPECIFICATION:

Model	Centrino Filtration Pump
Max. liquid flow rate	260ML Per min
Power	DC 24 V, Built-in Lithium battery
Motor	DC 24 V Pump
Hose barb	¼ Inch (6 mm)
Material contacts with	SS316L, Silicon Seal, PP Filtrate
Dimension	360mm(L)*160mm(W)*150mm(H)
Weight	3.7kg
Number of Branch	Triple Head
Filter holder	MF47 Funnel, PSU Funnel